<dummy: use COVER here>

USER INSTRUCTIONS—MANUEL DE L'UTILISATEUR—INSTRUCCIONES PARA EL USO

61301603331 Rev.A (12/15/2000)

Related Publications—Autres Publications—Publicaciones relacionadas

BIOSIGN Biological Indicator Instructions	350155
Installation & Operator Instructions for BIOSIGN Incubators	350156

IMPORTANT—IMPORTANT—IMPORTANTE

Special safety signs and instructions are used in this manual. Read them carefully before using the BIOSIGN EZ-VU TEST PACK. It is the responsibility of the user to follow all instructions.

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The following symbols with related notes appear in this manual.

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"Warning" notes alert the user to the possibility of personal injury.

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"Caution" notes alert the user to the possibility of damage to the equipment.

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"Notes" alert the user to pertinent facts and conditions.

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This manual contains proprietary information of Getinge/Castle, Inc., Rochester, New York USA.

It shall not be reproduced in whole or in part without the written permission of Getinge/Castle, Inc.

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
BIOSIGN EZ-VU TEST PACK

GENERAL DESCRIPTION

The BIOSIGN EZ-VU TEST PACK is a disposable biological indicator test pack. It is designed to provide both the challenge to steam sterilization processes that is presented by the AAMI 16-Towel Biological Indicator Challenge Test Pack, and the challenge to ethylene oxide (EO) sterilization processes that is presented by the AAMI Routine Test Pack for General-Purpose EO Sterilizers.

Each BIOSIGN EZ-VU TEST PACK contains a BIOSIGN Biological Indicator (BI) and a chemical indicator/recordkeeping strip. BIOSIGN BIs are vials that contain ampuled growth medium and filter paper with spores of *Bacillus stearothermophilus* and *Bacillus subtilis* var. *niger*. The ampuled growth medium has a pH indicator that changes color from RED to YELLOW if there is test organism growth. Two indicators show that a Test Pack has been exposed to sterilant. These indicators are on (1) the chemical indicator/recordkeeping strip, and (2) the label of the BIOSIGN BI. When the Test Pack is exposed to steam, a color appears next to STM. When the Test Pack is exposed to ethylene oxide, a color appears next to EO.

Separate BIOSIGN BIs are provided with each carton of BIOSIGN EZ-VU TEST PACKS. These are to be used as positive controls (BIs not processed in a sterilizer—see step 5 of BIOSIGN BIOLOGICAL INDICATOR (BI) INSTRUCTIONS below).

STORAGE

BIOSIGN EZ-VU TEST PACKS must be stored in a refrigerator or controlled room temperature ranging from 2° to 24°C (36° to 75°F). Store test packs away from excess heat, moisture, sterilants, or chemical vapors.

FREQUENCY OF USE

BIOSIGN EZ-VU TEST PACKS should be used in normal sterilization cycles:

- in Steam Sterilization—at least weekly. Daily use is recommended.
- in Ethylene Oxide (EO) Sterilization—in each sterilizer cycle.

Refer to the appropriate Association for the Advancement of Medical Instrumentation (AAMI) guidelines for the proper use of biological indicator test packs in health care facilities.

HOW TO USE THESE INSTRUCTIONS

- To use with steam sterilization, follow the steps under STEAM STERILIZATION.
- To use with ethylene oxide (EO) sterilization, follow the steps under ETHYLENE OXIDE (EO) STERILIZATION.
- To process the Indicator AFTER USE for either challenge, follow the steps under BIOSIGN BIOLOGICAL INDICATOR (BI) INSTRUCTIONS.

STEAM STERILIZATION

BIOSIGN EZ-VU TEST PACKS are not recommended for monitoring flash sterilization cycles at 132°C (270°F) or greater, also called “unwrapped method” cycles.

1. Load the sterilizer.


Before placing the BIOSIGN EZ-VU TEST PACK in a sterilizer, observe the ampule through the window. The ampule should contain red growth medium. If the ampule is dry or discolored, do not use this BIOSIGN EZ-VU TEST PACK for sterilizer monitoring. Dispose of it per institutional practice. Notify Getinge/Castle for a replacement.

Do not set any goods on top of the BIOSIGN EZ-VU TEST PACK. Otherwise, the test may be invalid.

2. Place one BIOSIGN EZ-VU TEST PACK window side up on the lowest shelf and above the drain.

3. Process the load. Use a normal cycle.

BURN HAZARD: *Wear heat-resistant gloves when removing the BIOSIGN EZ-VU TEST PACK. Materials are hot enough to burn unprotected skin.*

4. After the cycle is complete, remove the BIOSIGN EZ-VU TEST PACK. Let it cool to room temperature (at least 15 minutes).

Chemical indicators can only confirm that the load was processed. They cannot determine if the processing was enough to achieve sterility.

5. Open the end flaps of the BIOSIGN EZ-VU TEST PACK. Remove the BIOSIGN Biological Indicator (BI) and the chemical indicator/recordkeeping strip.
6. Check the chemical Indicator/recordkeeping strip. A color next to STM confirms that the sterilant has penetrated the pack. If the color did not appear, have the sterilizer checked for proper operation. Then repeat steps 1 through 6 using a new BIOSIGN EZ-VU TEST PACK.
7. Mark the appropriate information on the chemical indicator/recordkeeping strip. File the strip with the other sterilization processing records.
8. Let the BIOSIGN BI cool to room temperature. Then record the sterilizer load number and date on its label. A color next to STM confirms that it has been exposed to steam.

Biological indicator (BI) may burst if crushed when hot. Do not crush the ampule until the BI is cool. Let the BI cool for 15 minutes outside the sterilizer. Failure to cool may result in personal injury.


9. Activate the BI and incubate it. Follow the steps below (see BIOSIGN BIOLOGICAL INDICATOR (BI) INSTRUCTIONS).

ETHYLENE OXIDE (EO) STERILIZATION

EO EXPOSURE HAZARD: *Wear butyl gloves and a full face shield when exposure to EO is possible. Follow your institution's policies regarding EO safety.*

1. Load the sterilizer.
2. Remove the marked tab from each end flap as shown.

3. Set the BIOSIGN EZ-VU TEST PACK in the center of the load.

In small sterilizers, where load configuration does not permit placement of the BIOSIGN EZ-VU TEST PACK in the center of the load, the front portion of the sterilizer may be used.

4. Process the load. Use a normal cycle.

Wait for full aeration of loads processed by ethylene oxide before retrieving the BIOSIGN EZ-VU TEST PACK.

5. After aeration is complete, remove the BIOSIGN EZ-VU TEST PACK from the load.

Chemical indicators can only confirm that the load was processed. They cannot determine if the processing was enough to achieve sterility.


6. Open the end flaps of the BIOSIGN EZ-VU TEST PACK. Remove the BIOSIGN Biological Indicator (BI) and the chemical indicator/recordkeeping strip.
7. Check the chemical indicator/recordkeeping strip. A color next to EO confirms that it has been exposed to EO. If the color did not appear, have the sterilizer checked for proper operation. Then repeat steps 1 through 7 using a new BIOSIGN EZ-VU TEST PACK.
8. Mark the appropriate information on the chemical indicator/recordkeeping strip. File the strip with the other sterilization processing records.
9. Let the BIOSIGN BI cool to room temperature. Then record the sterilizer load number and date on its label. A color next to EO confirms that it has been exposed to EO.
10. Activate the BI and incubate it. Follow the steps below (see BIOSIGN BIOLOGICAL INDICATOR (BI) INSTRUCTIONS).

BIOSIGN BIOLOGICAL INDICATOR (BI) INSTRUCTIONS


Biological indicator (BI) may burst if crushed when hot. Do not crush the ampule until the BI is cool. Let the BI cool for 15 minutes outside the sterilizer. Failure to cool may result in personal injury.

1. Let the BI(s) cool at least 15 minutes, until the ampule cools to room temperature. Then snap the cap closed to seal uniformly against the collar as shown.

***RISK OF CUTS:** Close the cap first. Then use the ampule crusher. Using any other means or sequence to crush the ampule could result in personal injury.*

2. Activate the BI as shown. Use the ampule crusher built into the BIOSIGN incubator.

Failure to use the correct temperature incubator invalidates the testing and could lead to a false indication of sterility.

3. Place the activated BIs in the appropriate BIOSIGN incubator: 55°C for steam sterilization monitoring, or 37°C for ethylene oxide monitoring.

Continued incubation of positive BIs at 37°C for more than three days beyond the initial RED to YELLOW color change may result in reversion back to the original RED COLOR. This can be misunderstood and lead to a false interpretation of the test.

4. Growth is exhibited by a color change in the medium from RED to YELLOW. Check for growth daily until the recommended incubation time is achieved: 2 days (48 hours) at 55°C for steam sterilization monitoring, or 4 days (96 hours) at 37°C for ethylene oxide monitoring. Record the results in a record book. React to positive growth as soon as it is detected. Follow the policies of your organization.
5. On each day that you test a sterilizer, run a positive control (a BI not processed in a sterilizer) for each incubation temperature used. Follow the procedures outlined above for closing the cap, crushing the ampule, and incubating. Failure of the positive control to indicate growth is a serious condition that must be immediately addressed. Possible causes for non-growth include failure to activate the BI, a malfunctioning incubator, or inappropriate storage of BIs. Determine the cause for non-growth of positive controls and correct it.
6. Before disposal, steam sterilize all BIs for 30 minutes at 121°C (250°F), or incinerate them.

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© Castle and BIOSIGN are registered trademarks.—Castle et BIOSIGN sont des marques déposées.—
EZ-VU™ is a trademark of Getinge/Castle, Inc.—EZ-VU™ est une marque de fabrique de Getinge/Castle, Inc.—
Copyright © 2000 by Getinge/Castle, Inc.—Droits réservés © 2000 par Getinge/Castle, Inc.—
Copyright © 2000 por Getinge/Castle, Inc.